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## Preparation of electrocompetent cells

### OUTLINE:

Pre-culture (3 mL) → Culture (add the pre-cultured bacteria to the fresh medium at a ratio of 1:100 (e.g., 1 mL to 100 mL). Allow the culture to grow until O.D. 600 reaches 0.5–0.7 → Replace the culture medium with a glycerol solution.

**Procedure:** \*Perform all subsequent procedure within a tissue culture hood, whenever possible.

### Day0: Plating *E.coli* (evening)

- 1) Retrieve tubes containing competent *E.coli* cells (10 µL) from –80 °C (Deep freezer), and streak the cells on a LB plate without antibiotics.
- 2) Culture overnight at 37 °C.

### Day1\_morning: Advance preparation (NOT always necessary)

- 1) Autoclave 10–12 1.5 mL Eppendorf tubes. \*If you don't have autoclaved tubes, it is OK to use non-autoclaved ones.
- 2) Prepare 100 mL of LB in an Erlenmeyer flask (use larger 200–300 mL flask) and autoclave it. \*If you have an empty autoclaved flask, you can add the autoclaved LB to it when you are ready to culture.

### Day1\_evening: Pre-culture

- 3) Prepare approx. 3 mL of LB medium to a Spitz tube and inoculate the *E. coli* from a single colony in the medium.
- 4) Incubate overnight at 37 °C with shaking.

### Day2\_morning: Culture

- 1) Inoculate 1 mL of a fresh overnight *E.coli* culture into 100 mL of fresh LB medium (without antibiotics). \*The medium becomes just slightly cloudy.
- 2) Grow the cells at 37 °C with shaking to an OD600 of 0.5–0.7. \*The best results are obtained with cells that are harvested at early to mid-log phase (corresponding to a cell density of 4–5 x 10<sup>7</sup> cells/mL).

### Day2\_afternoon: Making competent cells

- 3) Chill the cells on ice for ~20 min. For all subsequent steps, keep the cells as close to 0°C as possible (in an ice/water bath) and chill all containers in ice before adding cells.
- 4) Transfer the cells to a clean, cold 50 mL conical tube and centrifuge at 4000 x g for 15 minutes at 4°C.

- 5) Carefully pour off and discard the supernatant. \*It is better to sacrifice yield by pouring off a few cells than to leave any supernatant behind. Gently resuspend the pellet in 100 mL of ice-cold 10% glycerol. Centrifuge at 4000 x g for 15 minutes at 4°C.
- 6) Carefully pour off and discard the supernatant. Resuspend the pellet in 50 mL of ice-cold 10% glycerol. Centrifuge at 4000 x g for 15 minutes at 4°C.
- 7) Carefully pour off and discard the supernatant. Resuspend the pellet in ~5 mL of ice-cold 10% glycerol. Centrifuge at 4000 x g for 15 minutes at 4°C.
- 8) Carefully pour off and discard the supernatant. Resuspend the cell pellet in a final volume of 0.4 mL of ice-cold 10% glycerol. Aliquot 100 µL of the suspension per 1.5 mL tube and keep the cells in a deep freezer (< -80 °C). For immediate use, proceed directly to electroporation.